Silicone-based soft denture liner
Weiche Silikonunterfütterung Für Prothesen
Materià de revèstiment souple en silicone pour prothèses dentaires

Manufacturer
Tokuyama Dental Corp.
38-9, Taïtuo 1-chome, Taïtuo-ku, Tokyo, Japan

EC Rep.
Tokuyama Europe GmbH
Ostrasse 10, 40211 Düsseldorf, Germany
Tel: +49-211-1754480 Fax: +49-211-357379

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North/South America/Middle East
J. MORITA USA, Inc.

9 Mason, Irvine, CA 92618, U. S. A.
Tel: (949) 581 9600 Fax: (949) 465 1095
Tel: (800) 831 3222 Email: info@jmoritausa.com

ENGLISH

SOFRELINER


Read all information, precautions, notes and materials before using the product.

PRODUCT DESCRIPTION AND GENERAL INFORMATION

SOFRELINER is an addition-cured silicone chaiside resilient lining material for dentures. Sofeliner has no odor or taste. It is provided in a new automatic cartridge delivery system, together with an adhesive primer accessories. The material cures in 5 minutes intraorally or in 20 minutes at room temperature (25°C/77°F).

CONTENTS

1. One cartridge of paste containing Base (27g) and Catalyst (27g), to be used with the TOKUYAMA DISPENSER II (Sold separately).....x1
 2. Primer (10 mL).....x2
 3. Brushes for Primer.....x2
 4. Spatula for spreading SOFRELINER after it has been applied to the denture.....x1
 5. TD II Mixing tip (S) or (XS) for use with paste cartridge and the TOKUYAMA DISPENSER II.....x1
 6. Steel Bur T 101.....x1
 7. Finishing Wheel.....x1
 8. Mandrel T.....x1
- * Note: Only the TOKUYAMA DISPENSER II can be used with SOFRELINER.

Refill

1. One cartridge of paste containing Base (27g) and Catalyst (27g)
2. Primer (10 mL)
3. TD II Mixing tip (S) (50 tips)
4. TD II Mixing tip XS (50 tips)
5. SOFRELINER FINISHING SYSTEM
 - * Steel Bur T 101.....x1
 - * Finishing Wheel.....x1
 - * Mandrel T.....x1

CLINICAL APPLICATIONS FOR SOFRELINER

1. Firm anhedged by sharp ridges of the alveolar bone and atrophy of the mucosa
2. Poor retention caused by extreme ridge atrophy
3. Major undercuts or maxillary/mandibular tori
4. Retention of overdentures before the final attachment is placed

STORAGE AND DISPOSAL

Keep SOFRELINER in a cool, dark place (0-25°C / 32-77°F) before using. After using SOFRELINER, the mixing tip must be replaced with the cartridge cap, which is attached to the cartridge at the time of purchase. DO NOT use the material if the expiration date has passed. To dispose the paste, dispense unused paste out from the cartridge and harden it before disposal. To dispose the Primer, remove the liquid from the bottle and soak the liquid in tissue paper before disposal.

CLINICAL PRECAUTIONS

1. SOFRELINER is intended for use only by licensed dental care professionals.
2. Use SOFRELINER only as instructed herein. DO NOT use SOFRELINER for any other purpose.
3. DO NOT use this product if the safety seals are broken or if the contents have been tampered with.
4. Primer contains methylene chloride, which is highly volatile. Avoid long-term exposure, ventilate the room during use, and close the cap tightly after use.
5. SOFRELINER is a silicone-based product. DO NOT use SOFRELINER on patients with a history of sensitivity to silicone. Always determine whether your patient is allergic to, or hypersensitive to silicone-based products. If a patient is uncertain whether or not he/she is sensitive to silicone based products, such as SOFRELINER, always perform a sensitivity test on the patient before using this product. Allow a small amount of SOFRELINER paste to set on a dispensing disc. Use this material as a test patch by placing the material with a sticking pliers on the patient's arm for at least 24 hours, monitoring the test patch area for sensitivity. Advise all patients of the possibility of sensitivity to silicone, and discuss the signs of silicone sensitivity with them. If a patient shows any signs of sensitivity, patient should be instructed to immediately contact their dentist for removal of SOFRELINER from the denture.
6. The use of examination gloves (plastic, vinyl) is recommended at all times to avoid the possibility of allergic reactions from SOFRELINER. The use of latex gloves may inhibit the cure of SOFRELINER. The SOFRELINER Primer comes in contact with the examination gloves, immediately dispose of the examination gloves and wash hands thoroughly with water.
7. AVOID CONTACT WITH EYES. In case of contact with eyes, thoroughly flush eyes with water, and immediately contact a physician.
8. Avoid contact of uncured chemical with skin. In case of contact with skin, clean effected area with rubbing alcohol, then thoroughly flush with water.
9. Avoid contact of unpolymerized materials with clothing. In case of contact with clothing, do not wipe the material. Allow it to set and peel it off.
10. DO NOT allow excess material to flow down the patient's trachea or esophagus. Aspirated or ingested SOFRELINER paste can cause asphyxiation or other injury, even though SOFRELINER paste is non-toxic. Immediately remove any excess lining material from the palatal area to avoid flow of material down the patient's throat. SOFRELINER is not radiopaque, and X-rays will not detect aspirated or ingested material.
11. Do not use SOFRELINER with any materials not listed in these instructions.
12. Clean the tip of the cartridge with a clean cloth or tissue each time before using a new mixing tip or replacing the cartridge cap.
13. The material cures in 5 minutes intraorally or in 20 minutes at room temperature (25°C / 77°F). However, SOFRELINER will not be fully cured until one day after the lining. Therefore, the softness of the material may change slightly during this period. The finishing procedure can be performed immediately after the extraction of the newly lined denture from the mouth or a gypsum model. The newly lined denture can be worn after the finishing procedure.
14. Always clean the spatula, brushes and cups with alcohol after each use.
15. Instruct the patients in accordance with the PATIENT PRECAUTIONS found in the last section of these instructions.

CLINICAL APPLICATION OF SOFRELINER

INTRAORAL LINING METHOD

To achieve the maximum effect using SOFRELINER the optimum thickness of the lining material should be 2-3 millimeters. The set of the material will not be affected if SOFRELINER is applied thicker than 3 millimeters in some areas. Areas which will require lining thicker than 3 millimeters should first be lined with TOKUSO REBASE. Tokuyama's chaiside hard denture lining material (sold separately), to obtain the proper thickness for SOFRELINER.

1. Remove the denture from the patient's mouth and clean the area to be lined with soap and water. Then air-dry the surface of the denture thoroughly.
2. Grind the tissue side of the denture to the required thickness with a bur. Bevel the borders of the denture to allow at least 2 millimeters of material to roll on the borders. Avoid setting the marginal line on the denture border.
3. After blowing the ground surface with strong air stream, pour the Primer into the plastic cup included in the SOFRELINER kit. Using the brush included in the SOFRELINER kit, apply the Primer evenly in a thin coat on the surface of the denture where SOFRELINER is to be lined. Do not miss any spot where the Primer should be applied, since this may cause poor adhesion.
4. Note: Attach the Drip-cap included in the SOFRELINER kit onto the rim of the Primer bottle to facilitate pouring.
- 5) Evenly dispense SOFRELINER paste on the tissue side of the denture. Spread the paste evenly with a spatula across all areas where the Primer was applied.
- 6) Suction any excess saliva from the patient's mouth 1 - 1 1/2 minutes after dispensing SOFRELINER onto the denture. Place the denture in the patient's mouth, and have the patient occlude very gently into centric occlusion to the proper vertical dimension. Check for any excess SOFRELINER in the patient's mouth, and remove immediately if necessary. The softness of the material will not be affected if SOFRELINER is applied thicker than 3 millimeters in some areas. Areas which will require lining thicker than 3 millimeters should first be lined with TOKUSO REBASE. Tokuyama's chaiside hard denture lining material (sold separately), to obtain the proper thickness for SOFRELINER.

1. Remove the denture from the patient's mouth and clean the area to be lined with soap and water. Then air-dry the surface of the denture thoroughly.
2. Grind the tissue side of the denture to the required thickness with a bur. Bevel the borders of the denture to allow at least 2 millimeters of material to roll on the borders. Avoid setting the marginal line on the denture border.
- 3) After blowing the ground surface with strong air stream, pour the Primer into the plastic cup included in the SOFRELINER kit. Using the brush included in the SOFRELINER kit, apply the Primer evenly in a thin coat on the surface of the denture where SOFRELINER is to be lined. Do not miss any spot where the Primer should be applied, since this may cause poor adhesion.
- 4) Note: Attach the Drip-cap included in the SOFRELINER kit onto the rim of the Primer bottle to facilitate pouring.
- 5) Evenly dispense SOFRELINER paste on the tissue side of the denture. Spread the paste evenly with a spatula across all areas where the Primer was applied.
- 6) Suction any excess saliva from the patient's mouth 1 - 1 1/2 minutes after dispensing SOFRELINER onto the denture. Place the denture in the patient's mouth, and have the patient occlude very gently into centric occlusion to the proper vertical dimension. Check for any excess SOFRELINER in the patient's mouth, and remove immediately if necessary. The softness of the material will not be affected if SOFRELINER is applied thicker than 3 millimeters in some areas. Areas which will require lining thicker than 3 millimeters should first be lined with TOKUSO REBASE. Tokuyama's chaiside hard denture lining material (sold separately), to obtain the proper thickness for SOFRELINER.

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